



## 2°C TO 8°C SMART VACCINE REFRIGERATOR BJY-3301

# 2°C TO 8°C PHARMACY REFRIGERATOR BJY-3301

SMART MEMORY MEDICAL REAGENT REFRIGERATOR



## SPECIFICATIONS

|                   |                                           |
|-------------------|-------------------------------------------|
| Model             | BJY-3301                                  |
| Capacity          | 360L                                      |
| Temperature Range | 2-8°C                                     |
| Outer size(W*D*H) | 680*660*1995mm                            |
| Inner size(W*D*H) | 1100*660*2065mm                           |
| Alt Name          | Smart Memory Medical Reagent Refrigerator |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)