



2°C TO 8°C SMART VACCINE REFRIGERATOR BJY-3302

2°C TO 8°C PHARMACY REFRIGERATOR BJY-3302

SMART MEMORY MEDICAL REAGENT REFRIGERATOR



SPECIFICATIONS

Model	BJY-3302
Capacity	660L
Temperature Range	2-8°C
Outer size(W*D*H)	580*490*1400mm
Inner size(W*D*H)	985*525*1375mm
Alt Name	Smart Memory Medical Reagent Refrigerator



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com