



-86°C ULTRA LOW TEMPERATURE CHEST FREEZER BJY-3401

-86°C ULTRA LOW TEMPERATURE FREEZER BJY-3401

CRYOGENIC TEST TUBE FREEZER



SPECIFICATIONS

Model	BJY-3401
Capacity	10L
Temperature Range	-86°C
Outer size(W*D*H)	530*525*855mm
Inner size(W*D*H)	190*290*190mm
Alt Name	Cryogenic test tube Freezer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com