



-145°C ULTRA LOW TEMPERATURE CHEST FREEZER BJY-3404

-150°C CRYOGENIC FREEZER BJY-3404

CRYOGENIC TEST TUBE FREEZER



SPECIFICATIONS

Model	BJY-3404
Capacity	10
Temperature Range	-145°C
Outer size(W*D*H)	560*575*1005mm
Inner size(W*D*H)	190*290*190mm
Alt Name	Cryogenic test tube Freezer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com