



PLATELET INCUBATOR BJY-3601

PLATELET INCUBATOR BJY-3601

BLOOD PLATELET SHAKING INCUBATOR



SPECIFICATIONS

Model	BJY-3601
Capacity	130L
Temperature Range	22±2
Outer size(W*D*H)	522*596*1365
Inner size(W*D*H)	432*500*843
Voltage	220V/50Hz
Cooling Type	Air cooling
Alt Name	Blood Platelet Shaking Incubator



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com