



**-60°C ULTRA LOW PORTABLE CAR FREEZER
BJY-403**

PORTABLE FREEZER BJY-403

PORTABLE CAR ULT FREEZER (AC)



High Energy Core
Famous compressor, more powerful cooling
Anti-vibration & Anti-shaking
Anti-vibration and anti-shake design for different road conditions
Dual-use
Support both DC and AC working,maximize the protection you want

SPECIFICATIONS

Model	BJY-403
Capacity	60L
Temperature Range	-60°C
Outer size(W*D*H)	816*488*575mm
Inner size(W*D*H)	650*320*390mm
Alt Name	Portable Car ULT Freezer (AC)

APPLICATIONS

For preservation of viruses, bacteria, cells, skin tissues, biological materials,etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com