



COMBINED REFRIGERATOR AND FREEZER BJJY-502

COMBINED REFRIGERATOR AND FREEZER BJY-502

VERTICAL DOUBLE TEMPERATURE REFRIGERATOR FREEZER



Multi-purpose

Refrigerator in 2-8°C & Freezer in -25°C, Partitioned Storage

Safe Keeping

There are two sets of independent refrigeration system for safety, efficiency and reliability

Environmental Protection

Refrigeration is 100% fluorine-free

Stable Performance

Famous compressor, more powerful cooling

SPECIFICATIONS

Model	BJY-502
Capacity	400L (280L + 120L)
Temperature Range	Refrigerator:2~8°C; Freezer:-25°C
Outer size(W*D*H)	770*715*1980mm
Inner size(W*D*H)	660*610*800/553*524*425mm
Alt Name	Vertical Double Temperature Refrigerator Freezer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com