



COMBINED REFRIGERATOR AND FREEZER BJJY-601

COMBINED REFRIGERATOR AND FREEZER BJY-601

VERTICAL DOUBLE TEMPERATURE REFRIGERATOR FREEZER



Multi-purpose
Refrigerator in 2-8°C & Freezer in -40°C, Partitioned Storage
Safe Keeping
There are two sets of independent refrigeration system for safety, efficiency and reliability
Environmental Protection
Refrigeration is 100% fluorine-free
Stable Performance
Famous compressor, more powerful cooling

SPECIFICATIONS

Model	BJY-601
Capacity	200 (120 + 80)
Temperature Range	Refrigerator:2~8°C; Freezer:-40°C
Outer size(W*D*H)	620*650*1680mm
Inner size(W*D*H)	530*490*540/425*435*409mm
Alt Name	Vertical Double Temperature Refrigerator Freezer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com