



COMBINED REFRIGERATOR AND FREEZER BJJY-701

COMBINED REFRIGERATOR AND FREEZER BJY-701

VERTICAL DOUBLE TEMPERATURE REFRIGERATOR FREEZER



Multi-purpose

Refrigerator at 2-8°C and freezer at -60°C, with partitioned storage.

Safe Keeping

Two independent refrigeration systems ensure safety, efficiency, and reliability.

Environmental Protection

The refrigeration system is 100% fluorine-free.

Stable Performance

Uses a well-known compressor for more powerful cooling.

SPECIFICATIONS

Model	BJY-701
Capacity	200 (120 + 80)
Temperature Range	Refrigerator: 2~8°C; Freezer: -60°C
Outer size(W*D*H)	620*650*1680mm
Inner size(W*D*H)	530*490*540/425*435*409mm
Alt Name	Vertical Double Temperature Refrigerator Freezer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com