



COMBINED REFRIGERATOR AND FREEZER BJJ-702

COMBINED REFRIGERATOR AND FREEZER BJY-702

VERTICAL DOUBLE TEMPERATURE REFRIGERATOR FREEZER



Multi-purpose

Refrigerator at 2-8°C and freezer at -60°C, with partitioned storage.

Safe Keeping

Two independent refrigeration systems ensure safety, efficiency, and reliability.

Environmental Protection

The refrigeration system is 100% fluorine-free.

Stable Performance

Uses a well-known compressor for more powerful cooling.

SPECIFICATIONS

Model	BJY-702
Capacity	400 (280 + 120)
Temperature Range	Refrigerator:2~8°C; Freezer:-60°C
Outer size(W*D*H)	770*715*1980mm
Inner size(W*D*H)	660*610*800/553*524*425mm
Alt Name	Vertical Double Temperature Refrigerator Freezer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com