



8°C TO 20°C PHARMACY REFRIGERATOR BGY-805

8°C TO 20°C PHARMACY REFRIGERATOR BJY-805

PHARMACY REFRIGERATOR 0 SERIES



Safety Design
Multiple alarm functions to protect your reagents in real time
Real Visibility
Built-in LED lighting for easy observation
Heating Glass Door
Ensure that the door body will not condense under high temperature and humidity environment
Constant Temperature
Forces air cooling system keeps the temperature stable and even and keeps the refrigerator dry

SPECIFICATIONS

Model	BJY-805
Capacity	660L
Temperature Range	8~20°C
Outer size(W*D*H)	1200*650*1895mm
Inner size(W*D*H)	1112*560*1232mm
Alt Name	Pharmacy refrigerator 0 series

APPLICATIONS

Used to store blood, reagents, medicines, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com