



## 8°C TO 20°C PHARMACY REFRIGERATOR BJY-806

# 8°C TO 20°C PHARMACY REFRIGERATOR BJY-806

## PHARMACY REFRIGERATOR 0 SERIES



### Safety Design

Multiple alarm functions to protect your reagents in real time

### Real Visibility

Built-in LED lighting for easy observation

### Heating Glass Door

Ensure that the door body will not condense under high temperature and humidity environment

### Constant Temperature

Forces air cooling system keeps the temperature stable and even and keeps the refrigerator dry

## SPECIFICATIONS

|                   |                                |
|-------------------|--------------------------------|
| Model             | BJY-806                        |
| Capacity          | 960L                           |
| Temperature Range | 8~20°C                         |
| Outer size(W*D*H) | 1200*745*1870mm                |
| Inner size(W*D*H) | 1080*610*1380mm                |
| Alt Name          | Pharmacy refrigerator 0 series |

## APPLICATIONS

Used to store blood, reagents, medicines, etc.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)