



MICRO CENTRIFUGE BJZ-1201

MICRO CENTRIFUGE BJZ-1201

MICROSCALE HIGH SPEED CENTRIFUGE



- 1.Unique airflow exchange technology ensures a temperature rise of only +10°C above room temperature,protecting temperature-sensitive samples.
2. 12x1.5/2.0ml centrifuge tubes, maximum relative centrifugal force: 15721xg (15000rpm)
- 3.Automatic opening of the centrifuge lid after centrifugation to prevent sample overheating and facilitate sample removal.
- 4.Equipped with an emergency lid-opening device for power outages, providing hidden safety and convenient operation.

SPECIFICATIONS

Model	BJZ-1201
Max. Speed	15000rpm
Max. RCF	15721xg
Max. Capacity	12x1.5/2ml
Total Power	≤ 50 W
Overall Machine Noise	≤ 56 dB(A)
Power Supply	AC 110-220 V, 50/60 Hz
Dimensions	23 x 24 x 12 cm
Net Weight	4.5 kg
Alt Name	Microscale High Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	name	Capacity (mL)	Max Speed (rpm)	Max RCF (xg)
1	Micro high-speed centrifuge	12 x 1.5/2	15,000	15,721
2	Fixed-angle rotor	12 x 1.5/2 (0.2 mL/0.5 mL adapters)	15,000	15,721

APPLICATIONS

Scientific Research, Clinical Testing, Biochemistry, Molecular Biology, Medical Labs, Universities, Research Institutions, Biotech Companies And All Sample Separation Scenarios.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com