



## BENCHTOP MICROPLATE CENTRIFUGE BJT-1301

# MICROPLATE CENTRIFUGE BJZ-1301

## MICROPLATE CENTRIFUGE



Inching centrifuge.  
Holds two 96-well microplates.  
Extremely small footprint, does not take up lab bench space.

## SPECIFICATIONS

|                        |                                        |
|------------------------|----------------------------------------|
| Model                  | BJZ-1301                               |
| Max. Speed             | 500-3000rpm                            |
| Max. RCF               | 600xg                                  |
| Max. Capacity          | 8x0.2ml                                |
| Sample throughput      | 2 x 96-hole PCR boards                 |
| Timing range           | Click/1s-99s/1min-99min                |
| Supported power supply | DC24V power adapter (AC85-265V50/60Hz) |
| Net weight             | 2.6kg                                  |
| Dimension              | 230*195*190(mm)                        |
| Alt Name               | Microplate Centrifuge                  |

## APPLICATIONS

Wide range of application: suitable for all kinds of standards with and without skirt edge

PCR microplate, 8x0.2mlpcr tube

Instant centrifugal function, instant acceleration to2500rpm, automatic brake stop time is fast

It can quickly remove the droplets hanging on the wall and is used before and after PCR experiments to significantly improve the experimental results

Maximum speed: 3000rpm, increment: 100rpm

Maximum relative centrifugal force: 600xg

Timing range: point move / 1-99 seconds / 1-99 minutes.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)