



BENCHTOP LOW SPEED CENTRIFUGE BJJZ-2201

LOW SPEED CENTRIFUGE BJZ-2201

BENCHTOP LOW SPEED CENTRIFUGE



- 1.Prevents safety risks from incorrect rotor installation.
- 2.Electronic lock + steel body + automatic unbalance stop.
- 3.Touch panel with precise speed/RCF/time settings, adjustable during operation.
- 4.Accommodates 0.2ml-750ml tubes, ELISA plates, blood bags, and culture flasks.
- 5.Timed, differential, and gradient centrifugation options.
- 6.Triple protection for virus and radioactive samples.
- 7.Noise level $\leq 60\text{dB (A)}$, suitable for long-term continuous work.

SPECIFICATIONS

Model	BJZ-2201
Max. Speed	6000rpm
Max. RCF	6321xg
Max. Capacity	4x500ml
Power Supply	AC 220 $\pm$ 22V 50Hz 10A
Net Weight	50KG
Total Power	800W
Noise Level	$\leq 60\text{dB(A)}$
Timer Range	1s-99min59s/1min-99h99min
Speed Accuracy	$\pm 10\text{r/min}$
Dimensions	540*470*350mm
Alt Name	Room Temperature Centrifuge

ACCESSORIES FOR PURCHASE

No	name	Rotor No.	Capacity (ml)	Max. Speed (rpm)	Max. RCF (*g)
1	Swing Out Rotor (Square Cup)	1	4x750ml (Square Cup)	4000	3580
2	Swing Out Rotor (Square Cup)	Adapter	4x500ml	4000	3580
3	Swing Out Rotor (Square Cup)	Adapter	4x250ml	4000	3580
4	Swing Out Rotor (Square Cup)	Adapter	4x3x100ml	4000	3580
5	Swing Out Rotor (Square Cup)	Adapter	4x7x50ml	4000	3580
6	Swing Out Rotor (Square Cup)	Adapter	4x14x5ml	4000	3580
7	Swing Out Rotor (Square Cup)	Adapter	4x14x10ml	4000	3580
8	Swing Out Rotor (Square Cup)	Adapter	4x24x5ml	4000	3580

9	Swing Out Rotor (Square Cup)	Adapter	4x30x1.5/2ml	4000	3580
10	Swing Out Rotor (Square Cup)	Microplate Rotor	4x4x96 hole	4000	3040

APPLICATIONS

This centrifuge is widely used in clinical medicine, biopharmaceuticals, cell culture, radioimmunoassay, and molecular laboratories. Optional biosafety accessories are available for handling hazardous infectious and radioactive samples.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com