



BENCHTOP LOW SPEED CENTRIFUGE BJT-2701

LOW SPEED CENTRIFUGE BJZ-2701

BENCHTOP LOW-SPEED CENTRIFUGE



SPECIFICATIONS

Model	BJZ-2701
Max. Speed	5000rpm
Max. RCF	3437xg
Max. Capacity	8x15ml+4x50ml
Power Supply	AC 220±22V 50Hz 10A
Net Weight	25KG
Total Power	500W
Noise Level	≤60dB(A)
Timer Range	1s-99min59s
Speed Accuracy	±10r/min
Dimensions	460*400*270mm
Alt Name	Benchtop low-speed centrifuge

ACCESSORIES FOR PURCHASE

No	name	Capacity	Max Speed	Max RCF
1	Fixed Angle Rotor	6 * 50ml	4000 rpm	2200xg
2	Swing Out Rotor	4 * 50ml + 8 * 15ml	4000 rpm	2200xg
3	Swing Out Rotor	24 * 5ml	4000 rpm	2200xg
4	Swing Out Rotor	8 * 15ml	4000 rpm	2200xg

APPLICATIONS

Public hospital lab, private clinic, blood bank: blood serum/plasma separation, urine routine test, PRP preparation, clinical biochemical sample separation.

Animal hospital & livestock testing lab: pet & livestock blood separation, animal serum analysis.

Biology/Medical college lab: gradient centrifugation, pulse spot centrifugation, biological reagent separation for academic research.

Medical ICL lab & food inspection lab: sample pre-treatment, suspension stratification, agricultural product testing.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com