



CONCENTRATOR FITTING-COLD TRAP BJZ-2801

COLD TRAP BATH BJZ-2801

CONCENTRATOR FITTING-COLD TRAP



SPECIFICATIONS

Model	BJZ-2801
Cold trap capacity	2L
Overall machine noise	≤20dB(A)
Cold trap type and size	200 mm x 120 mm
Dimension	250 x 430 x 360 (mm)
Cold trap temperature	-50°C
Power supply/power	220V/50Hz/0.4kW
Net weight	23kg
Alt Name	Concentrator fitting-Cold trap



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com