



PRP CENTRIFUGE BKA-1901

PRP CENTRIFUGE BKA-1901

PRP/ PORTABLE LOW SPEED CENTRIFUGE



1. Microprocessor control with Brushless Motor
2. Display the parameters of speed, timer, RCF, Error code information.
3. ABS Centrifuge body and light weight in 5KGS.
4. Safety Electrical lid interlock, machine won't run when lid open and error code appears.
5. Emergency lid lock release when power failure.
6. Window on lid for speed calibration and noise measurement.
7. Continuous operation, working continuously according to the set time.
8. Automatic troubleshooting with indication of error code and alarm on the screen.
9. Wide choice of rotors and adapters.
10. With CE and ISO9001:2008, ISO13485 certificate.

SPECIFICATIONS

Model	BKA-1901
Max. Speed	6000rpm
Max. RCF	2200xg
Max. Capacity	8x15ml
Noise(dBA)	≤58dB(A)
Time range	1-99minutes
Speed accuracy	±20r/min
Dimension(mm)(W*D*H)	310*270*210mm
Net weight(without rotor)	5KGS
Power supply	AC110~240V,50/60HZ,1A,0.04KW
Alt Name	PRP/ Portable Low Speed Centrifuge

ITEMS INCLUDED

No	Name	Max.Speed (rpm)	Max.Capacity	Max. RCF (xg)
1	Fixed Rotor	6000rpm	8x10ml	2100
2		6000rpm	8x15ml	2200
3		6000rpm	12x5ml(100mm blood tube)	2100



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com