



CELL WASHING CENTRIFUGE BKA-2001

CELL WASHING CENTRIFUGE BKA-2001

SERO CENTRIFUGE



1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.

SPECIFICATIONS

Model	BKA-2001
Max. Speed	4500rpm
Max. RCF	2000xg
Max. Capacity	12x7ml
Noise(dBA)	≤58dB(A)
Time range	1-99minutes
Speed accuracy	±20r/min
Dimension(mm)(W*D*H)	460*350*265mm
Net weight(without rotor)	23KGS
Power supply	AC110~240V,50/60HZ,1A,0.04KW
Alt Name	Sero Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max.Speed (rpm)	Max.Capacity	Max. RCF (xg)	Note
1	Swing rotor	4500	12x0.5/1.5/2.2ml	2000	To order
2	Swing rotor	3000	12x15ml	1000	To order

FEATURES

1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.
9. 10 kinds of Program(Memory) and 10 kinds of Acceleration&Deceleration for choice.If ask for more Programs please specially mention when place the order.
10. Automatic troubleshooting with indication of error code and alarm on the screen.
11. Wide choice of Rotors and adapters.
12. With CE and ISO9001:2008,ISO13485 certificate.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com