



PRP CENTRIFUGE BKA-2201

PRP CENTRIFUGE BKA-2201

PRP CENTRIFUGE



1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.

SPECIFICATIONS

Model	BKA-2201
Max. Speed	4000rpm
Max. RCF	3020xg
Max. Capacity	4x50ml
Noise(dBA)	≤58dB(A)
Time range	1-99minutes adjustable
Speed accuracy	±20r/min
Dimension(mm)(W*D*H)	460*585*360mm
Net weight(without rotor)	40KGS
Power supply	AC220-240V,50/60HZ,5A,0.5KW
Alt Name	PRP Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max.Speed (rpm)	Max.Capacity	Max. RCF (xg)	Note
1	Swing rotor	5000	4x50ml	4390	To order
2	Fixed rotor	4000	24x5ml(100mm Blood tube)	2600	To order

FEATURES

1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.
9. 10 kinds of Program(Memory) and 10 kinds of Acceleration&Deceleration for choice.If ask for more Programs please specially mention when place the order.
10. Automatic troubleshooting with indication of error code and alarm on the screen.
11. Wide choice of Rotors and adapters.
12. With CE and ISO9001:2008,ISO13485 certificate.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com