



MILK CENTRIFUGE BKA-2401

MILK CENTRIFUGE BKA-2401

GERBER MILK CENTRIFUGE



1. Microprocessor control with AC Brushless Motor
2. 5 Inches Colorful LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Heat release holes design reduce the temperature inside the centrifugal chamber on running.

SPECIFICATIONS

Model	BKA-2401
Max. Speed	2400 rpm
Max. RCF	1065xg
Max. Capacity	8PCS
Noise(dBA)	≤58dB(A)
Time range	1-99minutes
Temperature range(°C)	To 70°C
Speed accuracy	±20r/min
Dimension(mm)(W*D*H)	420*545*360mm
Net weight(without rotor)	43KGS
Power supply	AC220-240V,50/60HZ,10A,1KW
Alt Name	Gerber Milk Centrifuge

ACCESSORIES FOR PURCHASE

No	name	Capacity (ml)	Max. Speed (rpm)	Max. RCF (xg)
1	NO.1 Fixed rotor	8PCS	2400rpm	1065xg
2	NO.3 Fixed roto	4PCS	1500rpm	380xg

FEATURES

1. Microprocessor control with AC Brushless Motor
2. 5 Inches Colorful LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Heat release holes design reduce the temperature inside the centrifugal chamber on running.
8. Window on lid for speed calibration and noise measurement.
9. Continuous operation, working continuous according the set time.
10. 10 kinds of Program(Memory) and 10 kinds of Acceleration&Deceleration for choice.If ask for more Programs please specially mentioned when place the order.
11. Automatic troubleshooting with indication of error code and alarm on the screen.
12. Wide choice of Swing rotors and angle rotors and adapters.
13. With CE and ISO9001:2008,ISO13485 certificate.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com