



HEMATOCRIT CENTRIFUGE BKA-2501

HEMATOCRIT CENTRIFUGE BKA-2501

MICRO HEMATOCRIT AND MICRO TUBES CENTRIFUGE



1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,Program(Memory)information, Acceleration& 3.Deceleration information and Error code information.
4. ABS Centrifuge body and in light weight 8KGS.
5. Safety Electrical lid interlock,machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuous according the set time.
9. 10 kinds of Program(Memory) and 10 kinds of Acceleration&Deceleration for choice.If ask for more Programs please specially mentioned when place the order.
10. Automatic troubleshooting with indication of error code and alarm on the scre

SPECIFICATIONS

Model	BKA-2501
Max. Speed	16000rpm
Max. RCF	17800xg
Max. Capacity	24 Placers
Noise(dBA)	≤58dB(A)
Time range	1-99minutes
Speed accuracy	±20r/min
Dimension(mm)(W*D*H)	460*350*265mm
Net weight(without rotor)	8KGS
Power supply	AC110~240V,50/60HZ,3A,0.2KW
Alt Name	Micro Hematocrit and Micro Tubes Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max.Speed (rpm)	Max.Capacity	Max. RCF (xg)	Note
1	Hematocrit rotor	12000	24 Hematocrit rotor	14800	To order
2	Fixed rotor	16000	12*1.5/2.2ml	17800	To order

FEATURES

1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,Program(Memory)information, Acceleration& 3.Deceleration information and Error code information.
4. ABS Centrifuge body and in light weight 8KGS.
5. Safety Electrical lid interlock,machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuous according the set time.
9. 10 kinds of Program(Memory) and 10 kinds of Acceleration&Deceleration for choice.If ask for more Programs please specially mentioned when place the order.
10. Automatic troubleshooting with indication of error code and alarm on the screen.
11. 24 placers hematocrit rotor and micro rotor for choice.
12. With CE and ISO9001:2008,ISO13485 certificate



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com