



HEATING CENTRIFUGE BKA-2601

SPECIAL PURPOSE CENTRIFUGE BKA-2601

HEATING CENTRIFUGE



1. Microprocessor control with Brushless Motor
2. LED Display the parameters of speed ,timer,RCF,Acceleration&
3. Deceleration information and Error code information.
4. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
5. Gas hinge for lid drop protection.
6. Safety Electrical lid interlock,door can not open when rotor is running and The machine won't run when the lid opens and an error code appears.
7. Emergency lid lock release when power failure.
8. Window on lid for speed calibration and noise measurement.
9. Continuous operation, working continuously according to the set time.
10. 9 kinds of Acceleration&Deceleration for choice.

SPECIFICATIONS

Model	BKA-2601
Max. Speed	3500rpm
Max. RCF	3116xg
Max. Capacity	4x100ml
Noise(dBA)	58dB(A)
Time range	1-99minutes
Speed accuracy	±20r/min
Dimension(mm)(W*D*H)	540*670*480mm
Net weight(without rotor)	72KGS
Power supply	AC220~240V,50/60HZ,10A,1KW
Alt Name	Heating centrifuge

ITEMS INCLUDED

No	Name	Capacity (ml)	Max. Speed (rpm)	Max. RCF (xg)
1	NO.1 HEMATOCRIT ROTOR	4x100ml Adapter:4x50ml	3500rpm	3116 xg

FEATURES

1. Microprocessor control with Brushless Motor
2. LED Display the parameters of speed ,timer,RCF,Acceleration&
3. Deceleration information and Error code information.
4. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
5. Gas hinge for lid drop protection.
6. Safety Electrical lid interlock,door can not open when rotor is running and The machine won't run when the lid opens and an error code appears.
7. Emergency lid lock release when power failure.
8. Window on lid for speed calibration and noise measurement.
9. Continuous operation, working continuously according to the set time.
10. 9 kinds of Acceleration&Deceleration for choice.
11. Automatic troubleshooting with indication of error code and alarm on the screen.
12. 100ml crude oil separation rotor.

13. With CE and ISO9001:2015,ISO13485 certificate.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com