



FLOOR STANDING HIGH SPEED REFRIGERATED CENTRIFUGE BKB-1301

HIGH SPEED REFRIGERATED CENTRIFUGE BKB-1301

HIGH SPEED REFRIGERATED CENTRIFUGE



5-inch HD touch screen for super-fast parameter setting.

Multiple user-defined programs can be saved.

Up to 5 program groups can be selected to realize step centrifugation, diversify the experiment, and complete multiple centrifugation steps at one time.

11 acceleration and 12 deceleration rates. Among them, one ACC and DEC rate can be user-defined based on the experiment requirements.

The pre-cooling function can be started just by touching the button on the screen.

The running records and diagnostic records can be saved automatically, through which can effectively check the operation of the instrument and the centrifugation of each batch sample.

Low noise and maintenance-free brushless motor with long-service life.

SPECIFICATIONS

Model	BKB-1301
Max. Speed	10000rpm
Max. RCF	18300xg
Max. Capacity	6x1000 mL
Speed Accuracy	±10 rpm
Time Setting Range	1 min to 23 h 59 min
Temperature Setting Range	-20 °C to 40 °C
Temperature Accuracy	±2 °C
Compressor Unit	High Performance Compressor Unit; Environment-friendly refrigerant(R404a)
Noise	≤65 dB(A)
Power Supply	AC 220 V ± 22 V, 50 Hz
Power Consumption	5.5 kW
Dimensions (W x D x H)	710 x 820 x 1200 mm
Packing Dimensions (W x D x H)	920 x 1040 x 1520 mm
Net Weight	310kg
Alt Name	High Speed Refrigerated Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	10000	18300	6x250ml
2	Angle Rotor	8000	11620	6x500ml
3	Adapter	8000	11620	6x13x1.5ml
4				6x10x5ml
5				6x7x10ml
6				6x7x15ml(conical)
7				6x3x30ml
8				6x50ml(conical)
9				6x50ml
10				6x2x50ml(conical)
11				6x100ml

12				6x250ml
13	Angle Rotor	7000	11500	6x1000ml
14	Adapter	7000	11500	6x4x50ml
15				6x4x50ml(conical)
16				6x250ml
17				6x500ml
18	Swing Rotor Yoke	4000	4050	4x1000ml
19	Round buckets	4000	4050	4x1000ml
20	Adapter	4000	4050	4x500ml(conical)
21				4x500ml
22				4x300ml
23				4x250ml(conical)
24				4x1x200ml(blood bag)
25				4x3x100ml
26				4x5x50ml(conical)
27				4x5x50ml
28				4x14x15ml(conical)
29				4x14x15ml
30				4x37x7ml(RIA tube)
31				4x24x5/7ml(blood collection tube)
32				4x28x1.5ml
33	Swing Rotor Yoke	4000	4050	4x1000ml
34	Oval buckets	4000	4050	4x1000ml
35	Blood bag adapter	4000	4050	4x1000ml
36	Swing Rotor Yoke	4000	4050	4x1000ml
37	Microplate Buckets			4x1x96 well Deep 4x2x96 well



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com