



BENCHTOP HIGH SPEED CENTRIFUGE BKB-1701

HIGH SPEED CENTRIFUGE BKB-1701

HIGH SPEED CENTRIFUGE



Multiple user-defined programs can be saved.
9 acceleration and 10 deceleration rates, to achieve the best centrifugation result.
Equipped with multiple safety features, the high speed centrifuge machine includes over-speed protection, motor overheating safeguards, and a self-locking door lid to ensure the safety of both users and instruments.
Brushless DC motor, no pollution and maintenance-free, with long service life, high speed accuracy control.
Soft-touch lid closure for ergonomic lid locking.
Backlit LCD for visual monitoring parameters.
7075-T6 aviation grade aluminum alloy, hard oxidation, strong corrosion resistance.

SPECIFICATIONS

Model	BKB-1701
Max. Speed	16500rpm
Max. RCF	26054xg
Max. Capacity	6x50ml
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min
Noise	≤65dB(A)
Power Supply	AC220V±22V 50/60Hz
Power Consumption	400W
Dimension (WxDxH)	325x420x325(mm)
Packing Dimension (WxDxH)	385x640x330(mm)
Net Weight	20kg
Alt Name	High Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	16500	26054	24x1.5/2.2ml
2	Angle Rotor	15000	15845	12x5ml
3	Angle Rotor	13000	17320	12x10ml
4	Angle Rotor	11000	12840	6x50ml(round)
5	Angle Rotor	13500	18720	6x30ml
6	Angle Rotor	14000	22150	30x1.5/2.0ml
7	Angle Rotor	13000	16020	6x8x0.2mlPCR



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com