



BENCHTOP HIGH SPEED CENTRIFUGE BKB-1801

HIGH SPEED CENTRIFUGE BKB-1801

HIGH SPEED CENTRIFUGE



5 inch HD touch screen, simple interface, setting parameters, and running parameters can be displayed simultaneously.

Directly input the setting value for fast and convenient operation of the high speed centrifuge machine without adjusting the parameters back and forth many times. (compared to traditional keypad and knob control products)

Multiple user-defined programs can be saved.

Microprocessor control, DC brushless motor for stable and low noise operation.

Parameters can be set freely during a run

door lid interlock to improve the safety performance.

Compact with great performance.

10 Acceleration and 11 Deceleration rates, to achieve the best centrifugation result.

SPECIFICATIONS

Model	BKB-1801
Max. Speed	16500rpm
Max. RCF	18360xg
Max. Capacity	12x5ml
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min 59sec
Noise	≤65dB(A)
Power Supply	AC220V±22V 50Hz
Power Consumption	250W
Dimension (WxDxH)	320x370x235 (mm)
Packing Dimension (WxDxH)	530x410x290(mm)
Net Weight	11kg
Alt Name	High Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	16500	18360	12x1.5/2.2ml
2	Angle Rotor	16500	18360	18x0.5ml
3	Angle Rotor	12000	10142	12x5ml
4	Angle Rotor	13000	16170	24x1.5/2.2ml
5	Angle Rotor	15000	13800	3x8x0.2mlPCR Strip
6	Hematocrit rotor	12000	14800	12 capillaries
7	Hematocrit rotor	12000	14800	24 capillaries



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com