



BENCHTOP HIGH SPEED CENTRIFUGE BKB-1901

HIGH SPEED CENTRIFUGE BKB-1901

HIGH SPEED CENTRIFUGE



SPECIFICATIONS

Model	BKB-1901
Max. Speed	16000rpm
Max. RCF	21532xg
Max. Capacity	6x50ml
Speed Accuracy	±30/min
Time Setting Range	1min to 99min
Noise	≤62dB(A)
Power Supply	AC220V±22V 50/60Hz
Power Consumption	400W
Dimension (WxDxH)	330x390x300mm
Packing Dimension (WxDxH)	395x650x340mm
Weight	25kg
Alt Name	High Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	16000	17800	12x1.5/2.2ml
2	Angle Rotor	15000	15845	12x5ml
3	Angle Rotor	13000	17320	12x10ml
4	Angle Rotor	15000	21532	24x1.5/2.2ml
5	Angle Rotor	13500	18720	6x30ml
6	Angle Rotor	11000	12840	6x50ml



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com