



BENCHTOP LOW SPEED CENTRIFUGE BKB-2001

LOW SPEED CENTRIFUGE BKB-2001

LARGE CAPACITY LOW SPEED CENTRIFUGE



Touch interface for super-fast parameter setting compared with the traditional push-button operation.

RCF and Speed can be automatically switched for touch the key.

Equipped with a 7-inch full touch screen (which can be operated directly with gloves), the interface is simple and all parameters are clear at a glance.

Touch the button to automatically switch between RCF and Speed.

Two time setting ranges for optional: 1s-99min 59s and 1min to 99h59min

Three timing modes are available: start timing, speed timing, and continuous timing.

11 acceleration stages and 12 deceleration stages. And free parking time can be set.

Save up to 100 user-defined programs.

SPECIFICATIONS

Model	BKB-2001
Max. Speed	5500rpm
Max. RCF	4904xg
Max. Capacity	4x750ml
Speed Accuracy	±10r/min
Time Setting Range	1s to 99min59s or 1 min to 99h59min
Timing Mode	Start timing after spinning, Start timing after reach setting speed, Short spinning
Noise	≤65dB(A)
Power Supply	AC220V ±22V 50/60Hz
Power Consumption	900W
Dimension (WxDxH)	510x635x398mm
Packing Dimension (WxDxH)	660x770x450mm
Net Weight	76kg
Alt Name	Large Capacity Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	5500	4904	40x15ml
2	Swing Rotor Yoke	4000	3345	4x500ml
3	Round buckets			4x500ml
4	Swing Rotor Yoke	4000	3345	4x500ml
5	Round buckets(aerosol-tight)			
6	Adapter	4000	3345	4x1x250ml
7				4x2x100ml
8				4x4x50ml(conical/round)
9				4x7x20ml
10				4x9x15ml(conical/round)
11				4x11x15ml(conical)
12				4x14x10ml
13				4x19x5/7ml
14				4x36x1.5ml

15	Swing Rotor Yoke	4000	3345	4x500ml rectangular
16	Rectangular buckets			
17	Swing Rotor Yoke	4000	3345	4x500ml rectangular
18	Rectangular buckets(aerosol-tight)			
19	Adapter	4000	3345	4x1x250ml
20				4x3x100ml
21				4x6x50ml(conical/round)
22				4x12x15ml(conical)
23				4x12x15ml(round)
24				4x20x10ml
25				4x31x5/7ml
26				4x56x1.5ml
27	Swing Rotor Yoke	4000	3500	4x750ml
28	Round buckets			
29	Swing Rotor Yoke	4000	3500	4x750ml
30	Round buckets (aerosol-tight)			
31	Adapter	4000	3500	4x1x500ml(conical)
32		4000	3500	4x1x500ml(flat)
33		4000	3500	4x250ml(conical)
34		4000	3500	4x1x250ml(flat)
35		4000	3500	4x4x100ml
36		4000	3500	4x7x50ml(conical)
37		4000	3500	4x7x50ml
38		4000	3500	4x5x50ml(self-standing tube)
39		4000	3500	4x10x20ml
40		4000	3500	4x15x15ml(conical)
41		4000	3500	4x19x10ml(self-standing tube)
42		4000	3500	4x21x10ml
43		4000	3500	4x27x5/7ml
44		4000	3500	4x37x5ml(RIA tube)
45		4000	3500	4x56x1.5ml
46	Swing Rotor Yoke	4000	3500	4x750ml(Rectangular)
47	Rectangular buckets			
48	Adapter	4000	3500	4x9x50ml
49		4000	3500	4x9x50ml(conical)
50		4000	3500	4x19x15ml(conical)
51		4000	3500	4x19x15ml
52		4000	3500	4x30x10ml
53		4000	3500	4x38x5/7ml blood collection tube
54	Swing Rotor Yoke	4000	3500	4x2x96 well plate 4x4x96 well Plate total height≤67mm
55	Microplate Buckets			
56	Swing Rotor Yoke	4000	3370	Automatic Decapping
57	Hanging buckets	4000	3370	4x30pcs tube
58	Adapter	4000	3370	4x30pcs(blood collection tube)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com