



BENCHTOP LOW SPEED CENTRIFUGE BKB-2101

LOW SPEED CENTRIFUGE BKB-2101

LOW SPEED CENTRIFUGE



Various brackets and adapters are available.
Microprocessor control, DC brush-less motor.
RCF value can be calculated automatically.
Special damping structure to reduce vibration.
Countdown time displayed by seconds if the centrifugation time is less than 1min.

SPECIFICATIONS

Model	BKB-2101
Max. Speed	5000rpm
Max. RCF	4390xg
Max. Capacity	4x250ml
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min
Noise	≤65dB(A)
Power Supply	AC220V±22V 50/60Hz
Power Consumption	400W
Dimension (WxDxH)	450x545x380mm
Packing Dimension (WxDxH)	547x727x432mm
Weight	40kg
Alt Name	Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Swing Rotor Yoke	5000	4390	4x100ml
2	Brackets	5000	4390	4x100ml
3	Swing Rotor Yoke	4000	2810	
4	Brackets	4000	2810	8x50ml
5	Brackets			4x2x100ml
6	Brackets			32x15ml(10ml/5mlAdapter)
7	Microplate rotor	4000	2200	2x2x96 well
8	Microplate buckets	4000	2200	
9	Swing Rotor Yoke	4000	3000	4x250ml
10	Round buckets			4x250ml(flat)
11	Adapter	4000	3000	4x1x100ml
12				4x1x50ml(round)
13				4x1x50ml(conical)
14				4x4x20ml
15				4x4x15ml(conical)
16				4x6x15ml

17				4x7x10ml
18				4x10x5ml
19				4x8x1.5ml
20	Swing Rotor Yoke	4000	3000	4x250ml(Rectangular)
21	Swing rectangular buckets			
22	Adapter	4000	3000	4x250ml(flat)
23				4x4x50ml(conical)
24				4x7x15ml(conical)
25				4x9x10ml
26				4x9x5ml
27				4x15x1.5ml



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com