



BENCHTOP LOW SPEED CENTRIFUGE BKB-2201

LOW SPEED CENTRIFUGE BKB-2201

LOW SPEED CENTRIFUGE



Microprocessor control, DC brush-less motor.
Touch panel and LCD display.
Special damping structure to reduce vibration.
Integrated molding of housing.
Electronic door lock, centrifuge cannot work if the door is open and the door cannot be opened when it is working.
Countdown time displayed by seconds if the centrifugation time is less than 1min.

SPECIFICATIONS

Model	BKB-2201
Max. Speed	4000rpm
Max. RCF	2200xg
Max. Capacity	12x20ml
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min
Noise	≤ 62dB(A)
Power Supply	AC220V±22V 50/60Hz
Power Consumption	100W
Dimension (WxDxH)	456x330x266mm
Packing Dimension (WxDxH)	568x423x300mm
Weight	12kg
Alt Name	Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	4000	2200	12x20ml
2	Angle Rotor	4000	2200	18x10ml
3	Angle Rotor	4000	2000	4x50ml
4	Angle Rotor	4000	2200	24x10ml



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com