



BENCHTOP LOW SPEED CENTRIFUGE BKB-2301

LOW SPEED CENTRIFUGE BKB-2301

LOW SPEED CENTRIFUGE



Variety of brackets accommodate different capacity tubes.
RCF value can be calculated automatically.
Save up to 25 user-defined programs.
Microprocessor control and DC brushless motor.
Special damping structure can reduce vibration.
Close the centrifuge cover with a soft touch cover. If the door is open, the centrifuge cannot work, and the door cannot be opened while it is working.
Touch panel and LCD display.

SPECIFICATIONS

Model	BKB-2301
Max. Speed	4200rpm
Max. RCF	2760xg
Max. Capacity	12x20ml
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min
Noise	≤65dB(A)
Power Supply	AC220V±22V 50/60Hz
Power Consumption	250W
Dimension (WxDxH)	380x500x300mm
Packing Dimension (WxDxH)	465x660x350mm
Weight	19kg
Alt Name	Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	4000	2200	18x10ml
2	Angle Rotor	4000	2200	12x20ml
3	Swing rotor yoke(round disc)	4200	2760	4x50ml
4	Brackets	4200	2760	4x50ml(round)
5	Adapter	4200	2760	4x15ml(conical)
6	Swing rotor yoke(round disc)	4200	2760	12x15ml
7	Brackets	4200	2760	12x15ml(round)
8	Adapter	4200	2760	12x1.5ml
9				12x5ml
10				12x10ml
11	Brackets	4200	2760	12x2x5ml
12	Adapter	4200	2760	12x2x1.5ml
13	Adapter	4200	2760	12x2x0.5ml
14	Swing rotor yoke(round disc)	4200	2760	8x10ml
15	Brackets	4200	2760	8x10ml
16	Brackets	4200	2760	8x2x5ml



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com