



BENCHTOP LOW SPEED CENTRIFUGE BKB-2401

LOW SPEED CENTRIFUGE BKB-2401

LOW SPEED CENTRIFUGE



A variety of brackets accommodate different capacity tubes.
Microprocessor control and DC brushless motor.
Touch panel and LCD display.
RCF value can be calculated automatically.
Special damping structure to reduce vibration.
Centrifuge lid with soft-touch lid closure, centrifuge cannot work if the door is open and the door cannot be opened when it is working.
Save up to 25 user-defined programs.

SPECIFICATIONS

Model	BKB-2401
Max. Speed	5000rpm
Max. RCF	4030xg
Max. Capacity	16x15ml
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min
Noise	≤65dB(A)
Power Supply	AC220V±22V 50/60Hz
Power Consumption	250W
Dimension (WxDxH)	380x500x300mm
Packing Dimension (WxDxH)	465x660x350mm
Weight	19kg
Alt Name	Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	5000	4030	24x10ml
2	Adapter	5000	4030	24x1.5ml
3	Swing Rotor Yoke	4200	2760	4x50ml
4	Brackets	4200	2760	4x50ml (conical/round)
5	Brackets	4200	2760	4x4x1.5ml(conical/round) (10ml/5mlAdapter)
6	Brackets	4200	2760	4x6x5ml(blood collection tube)
7	Adapter	4200	2760	4x6x1.5ml
8	Adapter	4200	2760	4x6x0.5ml



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com