



BENCHTOP LOW SPEED CENTRIFUGE BKB-2501

LOW SPEED CENTRIFUGE BKB-2501

LOW SPEED CENTRIFUGE



Touch panel and LCD display.
RCF value can be calculated automatically.
Save up to 25 user-defined programs.
Microprocessor control and DC brushless motor.
Special damping structure to reduce vibration.
Centrifuge lid with soft-touch lid closure. A centrifuge cannot work if the door is open and the door cannot be opened when it is working.
A variety of brackets accommodate different capacity tubes.

SPECIFICATIONS

Model	BKB-2501
Max. Speed	6000rpm
Max. RCF	5120xg
Max. Capacity	16x15ml
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min
Noise	≤65dB(A)
Power Supply	AC220V±22V 50/60Hz
Power Consumption	250W
Dimension (WxDxH)	380x500x300mm
Packing Dimension (WxDxH)	465x660x350mm
Weight	19kg
Alt Name	Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	6000	5120	12x15ml
2	Swing Rotor Yoke	4200	2760	4x50ml
3	Brackets	4200	2760	4x50ml(conical/round)
4	Brackets			4x4x15ml(conical/round) (10ml/5mlAdapter)
5	Brackets			4x6x5ml(Blood collection tube)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com