



BENCHTOP LOW SPEED REFRIGERATED CENTRIFUGE BKB-2701

LOW SPEED REFRIGERATED CENTRIFUGE BKB-2701

LOW SPEED REFRIGERATED CENTRIFUGE



It adopts a double-layer structure. A compact footprint (400x470x350mm) saves valuable bench space.

With 5 inch HD touch screen control, it is easy-operating and intuitive. The setting parameters and running parameters at a glance.

Multiple adapters available, accommodate 2- 5ml vacuum blood collection tube, 8 - 10ml nucleic acid centrifuge tube as well as other 10ml tubes.

Multiple protection, include over-speed, over-temperature, lock and so on, to ensure user's safety.

SPECIFICATIONS

Model	BKB-2701
Max. Speed	5000rpm
Max. RCF	3354xg
Max. Capacity	8x10ml
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min
Temperature Setting Range	-20 to 40°C
Temperature Accuracy	±2.0°C
Compressor Unit	High Performance Compressor Unit; Environment-friendly Refrigerant(R134a)
Noise	≤60dB(A)
Power Supply	AC220V±22V 50/60Hz
Power Consumption	450W
Dimension (WxDxH)	330x450x330mm
Packing Dimension (WxDxH)	720x520x550mm
Weight	44kg
Alt Name	Low Speed Refrigerated Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Swing Rotor (whole set , Yoke & buckets)	4500	3080	swing rotor
2	Adapter	4500	3080	4x6x5/7ml
3	Adapter	4500	3080	4x4x10ml
4	Angle Rotor	5000	3555	12x15ml
5	Angle Rotor	5000	4030	24x10ml



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com