



BENCHTOP LOW SPEED REFRIGERATED CENTRIFUGE BKB-2801

LOW SPEED REFRIGERATED CENTRIFUGE BKB-2801

LOW SPEED REFRIGERATED CENTRIFUGE



With a unique sucker type foot, it can be sucked on the table tightly. With 5 inch HD touch screen control, it is easy-operating and intuitive. The setting parameters and running parameters at a glance. Multiple adapters available, it's a compact, versatile and indispensable general purpose centrifuge. An environmentally friendly compressor is adopted, which has a good refrigeration effect and no pollution. Multiple protection, include over-speed, over-temperature, lock and so on, to ensure user's safety.

SPECIFICATIONS

Model	BKB-2801
Max. Speed	5000rpm
Max. RCF	3555xg
Max. Capacity	4x4x15ml or nucleic acid tube
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min
Temperature Setting Range	-20 to +40°C
Temperature Accuracy	±1°C
Compressor Unit	High Performance Compressor Unit; Environment-friendly Refrigerant(R134a)
Noise	≤65dB(A)
Power Supply	AC220V±22V 50Hz
Power Consumption	500W
Dimension (WxDxH)	415x625x320mm
Packing Dimension (WxDxH)	960x530x520mm
Weight	55kg
Alt Name	Low Speed Refrigerated Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	5000	3555	12x15ml
2	Swing Rotor Yoke	4200	2958	
3	Brackets	4200	2958	4x50ml
4	Adapter	4200	2958	4x15ml(conical)
5	Brackets	4200	2958	4x4x15ml
6	Adapter	4200	2958	4x4x1.5ml
7				4x4x5ml
8	Brackets	4200	2958	4x6x5ml



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com