



FLOOR STANDING HIGH SPEED REFRIGERATED CENTRIFUGE BKB-3601

HIGH SPEED REFRIGERATED CENTRIFUGE BKB-3601

HIGH SPEED REFRIGERATED CENTRIFUGE



12 kinds rotors options, which can provide either high-speed or large-capacity centrifugal operation.
Over-speed, Over-temperature protection and door lid self-locked, stainless steel chamber for safe operation.
Microprocessor control, AC brush-less motor, without carbon dust pollution and service life is long.
9 Acceleration and 10 Deceleration rates (NO.0 is free parking), 3 levels damping, to achieve best centrifugal result.
The automatic rotor identification system against over-speed.
The environmentally friendly refrigerant can keep the temperature below -4°C at max. speed.
Touch panel, LCD display, programmable.

SPECIFICATIONS

Model	BKB-3601
Max. Speed	20500rpm
Max. RCF	29200xg
Max. Capacity	4x750ml
Speed Accuracy	±10r/min
Temperature Setting Range	-20°C to +40°C
Temperature Accuracy	±1°C
Time Setting Range	1min to 99min 59sec or HOLD
Compressor Unit	High Performance Compressor Unit; Environment-friendly Refrigerant(R404a)
Noise	≤65dB(A)
Power Supply	AC220V±22V 50Hz
Power Consumption	1500W
Dimension (WxDxH)	550x650x840mm
Packing Dimension (WxDxH)	756x856x1320mm
Weight	110kg
Alt Name	High Speed Refrigerated Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	20500	29200	12x1.5/2.2ml
2	Angle Rotor	16500	19175	12x5ml
3	Angle Rotor	15000	23120	12x10ml
4				6x30ml
5	Angle Rotor	13000	17940	6x50ml(round)
6	Angle Rotor	13000	17940	6x50ml(conical)
7	Swing Rotor Yoke	4000	3500	4x750ml
8	Round buckets			
9	Swing Rotor Yoke	4000	3500	4x750ml
10	Round buckets(aerosol-tight)			

11	Swing Rotor Yoke	4000	3345	4x1x96 well deep 4x2x96 well
12	Microplate buckets			
13	Swing Rotor Yoke	4000	3345	4x500ml
14	Round buckets			
15	Swing Rotor Yoke	4000	3345	4x500ml
16	Round buckets(aerosol-tight)			
17	Swing Rotor Yoke	4000	3345	4x500ml rectangular
18	Rectangular buckets			
19	Swing Rotor Yoke	4000	3345	4x500ml rectangular
20	Rectangular buckets(aerosol-tight)			
21	Angle Rotor	16500	26050	24x1.5/2.2ml
22	Angle Rotor	12000	15805	6x85ml
23	Angle Rotor	10000	11960	6x100ml
24	Angle Rotor	15000	24036	8x15ml(conical)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com