



MINI CENTRIFUGE BKB-4302

MINI CENTRIFUGE BKB-4302

MINI CENTRIFUGE



Voiceless technology, running with low noise
Equipped with a brushless motor, the service life is long.
Equipped with a 6 place 1.5/2mL tube rotor as well as a 2-place PCR 8-strip tube rotor

AC100 to 242V, global application.

Each Cence small centrifuge for sale is engineered to be durable and safe for users, providing reliability you can count on.

WTL mini centrifuge takes simple one-step operation to a new level. Simply place your sample in the centrifuge, close the lid, and the small centrifuge for sale will begin to rotate without the press of a button.

SPECIFICATIONS

Model	BKB-4302
Max. Speed	6000rpm
Max. RCF	2000xg
Max. Capacity	6x1.5ml/0.5ml/0.2ml
Noise	≤53dB
Power Supply	AC100 to 242V
Power Consumption	30W
Dimension (WxDxH)	160x200x123(mm)
Packing Dimension (WxDxH)	205x265x160(mm)
Weight	1.2Kg
Alt Name	Mini Centrifuge

ITEMS INCLUDED

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Rotor	4000	900	6x1.5ml/0.5ml/0.2ml 2x8x0.2ml rotor



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com