



DENTAL CGF CENTRIFUGE BKB-4601

PRP CENTRIFUGE BKB-4601

LOW SPEED CGF CENTRIFUGE FOR DENTISTRY



The microcomputer-controlled centrifuge, driven by a DC brushless motor, operates smoothly and quietly for an enhanced user experience.

The centrifuge allows parameter adjustments during operation, eliminating the need for shutdown and ensuring seamless performance.

When the countdown time is less than one minute, the centrifuge for PRP for sale is displayed in seconds.

The centrifuge is equipped with a fault self-diagnosis system, capable of detecting issues such as overspeed, speed loss, and door cover protection to enhance operational safety.

With its sleek and stylish appearance, the centrifuge is compact and incorporates special noise treatment, ensuring an overall noise level of under 62 dB(A).

The centrifuge features a 5-inch HD touchscreen control, offering an intuitive i

SPECIFICATIONS

Model	BKB-4601
Max. Speed	4000rpm
Max. RCF	1980xg
Max. Capacity	8x15ml
Speed Accuracy	±30r/min
Time Setting Range	1min~99min
Noise	≤62 dB(A)
Power Supply	AC220v±22V 50/60Hz
Power Consumption	100W
Dimension (WxDxH)	320x370x235mm
Packing Dimension (WxDxH)	530x410x290mm
Weight	11Kg
Alt Name	Low Speed CGF Centrifuge for Dentistry

ITEMS INCLUDED

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle rotor	4000	1980	8x15ml

FEATURES

The microcomputer-controlled centrifuge, driven by a DC brushless motor, operates smoothly and quietly for an enhanced user experience.

The centrifuge allows parameter adjustments during operation, eliminating the need for shutdown and ensuring seamless performance.

When the countdown time is less than one minute, the centrifuge for PRP for sale is displayed in seconds.

The centrifuge is equipped with a fault self-diagnosis system, capable of detecting issues such as overspeed, speed loss, and door cover protection to enhance operational safety.

With its sleek and stylish appearance, the centrifuge is compact and incorporates special noise treatment, ensuring an overall noise level of under 62 dB(A).

The centrifuge features a 5-inch HD touchscreen control, offering an intuitive interface where set and running parameters are conveniently displayed on the same screen.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com