



CYTO CENTRIFUGE BKB-4701

CYTO CENTRIFUGE BKB-4701

CYTOSPIN CENTRIFUGE



Microprocessor control and DC brushless motor
Touch panel and LCD display. RCF can be calculated automatically.
Save up to 25 self-defined programs. This cytocentrifuge features a separate short spin key for quick centrifugation.
Special damping structure to reduce vibration, ensuring consistent results in cell debris centrifugation.
Centrifuge lid with soft-touch lid closure. Centrifuge cannot work if the door is open and the door cannot be opened when it's working.
Several cell smear rotors available. It can smear various fluid to achieve the best centrifugation effect.

SPECIFICATIONS

Model	BKB-4701
Max. Speed	3000rpm
Max. RCF	890xg
Max. Capacity	12-place
Funnel Capacity	0.1ml-0.5ml
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min
Noise	≤60dB(A)
Power Supply	AC220V ±22V 50/60Hz
Power Consumption	260W
Dimension (WxDxH)	380x500x300mm
Packing Dimension (WxDxH)	465x660x350mm
Weight	19kg
Alt Name	Cytospin Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Swing buckets cell smear rotor	2000	545	6 buckets
2	Cell smear rotor	3000	730	6 positions
3	Cell smear rotor	3000	890	12 positions
4	Box type rotor	2000	405	6 positions



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com