



AUTOMATIC DECAPPING CENTRIFUGE BKB-4801

AUTOMATIC DECAPPING CENTRIFUGE BKB-4801

LOW SPEED AUTOMATED DECAPPING CENTRIFUGE



Microcomputer control, LCD display, Variable frequency AC brushless motor.
The rotor can be sterilized under high temperature and pressure.
The vacuum blood collection tube is automatically capped with more than 90%, convenient and quick. Can effectively prevent hemolysis caused by the vibration when manually removing the cap.
Multiple self-defined programs can be saved

SPECIFICATIONS

Model	BKB-4801
Max. Speed	4000rpm
Max. RCF	3030xg
Max. Capacity	80pcs(13x75/100mm)
Speed Accuracy	±10r/min
Time Setting Range	1sec to 99min59sec
Noise	≤65dB(A)
Power Supply	AC220V ±22V 50Hz
Power Consumption	700W
Dimension (WxDxH)	470x560x380mm
Packing Dimension (WxDxH)	547x727x422mm
Net Weight	50kg
Alt Name	Low Speed Automated Decapping Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Swing buckets cell smear rotor	4000	3030	80 pcs(13x75~100)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com