



AUTOMATIC DECAPPING CENTRIFUGE BKB-5101

AUTOMATIC DECAPPING CENTRIFUGE BKB-5101

LOW SPEED AUTOMATED DECAPPING CENTRIFUGE



LCD touch screen with simple & clear interface for viewing and setting parameters of Speed, Time, RCF, Rotors, Program (Memory), Acceleration/Deceleration rate. Observations window on lid for speed calibration and noise measurement. Microprocessor control with frequency conversion A.C. motor drive. 10 acceleration rates and 11 braking rates, to effectively avoid secondary sedimentation, with better centrifugation effect.

SPECIFICATIONS

Model	BKB-5101
Max. Speed	4000rpm
Max. RCF	3030xg
Max. Capacity	80pcs (13x75/100mm tubes)
Speed Accuracy	±10r/min
Time Setting Range	1sec to 99min59sec
Temperature Setting Range	-20 to 40°C
Temperature Accuracy	±1°C
Compressor Unit	High Performance Compressor Unit; Environment-friendly refrigerant(R404a)
Noise	≤62dB(A)
Power Supply	AC220V ±22V 50Hz
Power Consumption	1500W
Dimension (WxDxH)	710x560x450mm
Packing Dimension (WxDxH)	1022x712x630mm
Net Weight	95kg
Alt Name	Low Speed Automated Decapping Centrifuge

ITEMS INCLUDED

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Swing Rotor	4000	3030	80 pcs (13x75~100)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com