



AUTOMATIC DECAPPING REFRIGERATED CENTRIFUGE BKB-5201

AUTOMATIC DECAPPING CENTRIFUGE BKB-5201

LARGE CAPACITY REFRIGERATED AUTOMATED DECAPPING CENTRIFUGE



Compatible with 13x75/100mm vacuum blood collection tube, more than 90% tubes more be successfully decapped.
Automatic calculation of RCF value, electronic door lock, strong safety performance

SPECIFICATIONS

Model	BKB-5201
Max. Speed	4000rpm
Max. RCF	3370xg
Max. Capacity	120pcs (13x75/100mm tubes)
Speed Accuracy	±10r/min
Temperature Setting Range	15 to 25°C
Temperature Accuracy	±2°C
Compressor Unit	High Performance Compressor Unit; Environment-friendly refrigerant(R134a)
Time Setting Range	1sec to 99min59sec
Noise	≤62dB(A)
Power Supply	AC220V ±22V 50Hz
Power Consumption	1800W
Dimension (WxDxH)	560x660x850mm
Packing Dimension (WxDxH)	756x856x1320mm
Net Weight	130kg
Alt Name	Large Capacity Refrigerated Automated Decapping Centrifuge

ITEMS INCLUDED

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Swing Rotor	4000	3370xg	120pcs (13x75~100)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com