



## AUTOMATIC DECAPPING REFRIGERATED CENTRIFUGE BKB-5301

# AUTOMATIC DECAPPING CENTRIFUGE BKB-5301

LARGE CAPACITY REFRIGERATED AUTOMATED DECAPPING CENTRIFUGE



DC brushless motor drive, low noise, maintenance free, no dust.  
Applied to diameter 13mm, height 100mm and 75mm tubes. (also suitable for diameter 16mm vacuum blood collection tube, adapter need customized)

## SPECIFICATIONS

|                           |                                                                           |
|---------------------------|---------------------------------------------------------------------------|
| Model                     | BKB-5301                                                                  |
| Max. Speed                | 4000rpm                                                                   |
| Max. RCF                  | 3370xg                                                                    |
| Max. Capacity             | 120pcs (13x75/100mm tubes)                                                |
| Speed Accuracy            | ±10r/min                                                                  |
| Temperature Setting Range | 0 to 35°C                                                                 |
| Temperature Accuracy      | ±1°C                                                                      |
| Compressor Unit           | High Performance Compressor Unit; Environment-friendly refrigerant(R134a) |
| Time Setting Range        | 1sec to 99min59sec                                                        |
| Noise                     | ≤62dB(A)                                                                  |
| Power Supply              | AC220V ±22V 50Hz                                                          |
| Power Consumption         | 2000W                                                                     |
| Dimension (WxDxH)         | 565x670x850mm                                                             |
| Packing Dimension (WxDxH) | 756x856x1320mm                                                            |
| Net Weight                | 140kg                                                                     |
| Alt Name                  | Large Capacity Refrigerated Automated Decapping Centrifuge                |

## ITEMS INCLUDED

| No | Name             | Max. Speed | Max. RCF | Max. Capacity |
|----|------------------|------------|----------|---------------|
| 1  | Swing Rotor Yoke | 4000       | 3370     |               |
| 2  | Swing Bucket     | 4000       | 3370     |               |
| 3  | Adapter          | 4000       | 3370     | 4x30x7ml/5ml  |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)