



AUTOMATIC DECAPPING CENTRIFUGE BKB-5401

AUTOMATIC DECAPPING CENTRIFUGE BKB-5401

LARGE CAPACITY LOW SPEED AUTOMATED DECAPPING CENTRIFUGE



Touch Screen for easier operation. Driven by AC frequency-conversion motor, with functions of accurate speed control, fast acceleration and deceleration.
The decap rate not less than 90%
Avoid re-mixing of blood after separation, prevent bacterial infection

SPECIFICATIONS

Model	BKB-5401
Max. Speed	4000rpm
Max. RCF	3850xg
Max. Capacity	150pcs (13x75/100mm tubes)
Speed Accuracy	±10r/min
Time Setting Range	1sec to 99min59sec
Noise	≤65dB(A)
Power Supply	AC220V ±22V 50Hz
Power Consumption	2000W
Dimension (WxDxH)	700x880x930mm
Packing Dimension (WxDxH)	820x1000x1080mm
Net Weight	220kg
Alt Name	Large Capacity Low Speed Automated Decapping Centrifuge

ITEMS INCLUDED

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Swing Rotor Yoke	4000	3850	150 pcs (13x75/100)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com