



AUTOMATIC DECAPPING REFRIGERATED CENTRIFUGE BKB-5501

AUTOMATIC DECAPPING CENTRIFUGE BKB-5501

AUTOMATIC DE-CAPPING REFRIGERATED CENTRIFUGE



Automatic troubleshooting with indication of error code and alarm on the screen.
Special damping set releasing the vibration.
Multiple groups of user-defined programs can be saved.
Applied to blood collection tubes and more than 90% tubes can be automatically decapped.
Suitable for vacuum blood collection tubes (13x75/100mm)

SPECIFICATIONS

Model	BKB-5501
Max. Speed	4000rpm
Max. RCF	3850xg
Max. Capacity	150pcs (13x75/100mm)
Speed Accuracy	±10r/min
Temperature Setting Range	-20 to 40°C
Temperature Accuracy	±1°C
Compressor Unit	High Performance Compressor Unit; Environment-friendly refrigerant(R404a)
Time Setting Range	1sec to 99min59sec
Noise	≤65dB(A)
Power Supply	AC220V ±22V 50Hz
Power Consumption	3000W
Dimension (WxDxH)	700x880x930mm
Packing Dimension (WxDxH)	820x1000x1080mm
Net Weight	260kg
Alt Name	Automatic De-capping Refrigerated Centrifuge

ITEMS INCLUDED

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Swing Rotor Yoke	4000	3850	150 pcs (13x75/100)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com