



OIL TEST CENTRIFUGE BKB-5601

OIL TEST CENTRIFUGE BKB-5601

OIL TESTING CENTRIFUGE



LCD of centrifuge machine for crude oil testing to show rotor No., speed, time, temperature, RCF, Acc/Dec rate, and trouble information/error code.

A circular heater for fast heating and temperature uniformity realizes the functions of temperature controlling and temperature holding, which meet the temperature requirements for the oil industry, making it a perfect centrifuge machine for crude oil testing.

The frequency-conversion motor of centrifuge machine for crude oil testing with high torque is characterized by stable and fast acceleration/deceleration as well as long service life, maintenance-free, and low noise for stable operation.

The centrifuge for crude oil testing is equipped with a specially designed tube, making it resistant to damage during operation. The oil centrifuge is a widely utilized tool in Daqing Oilfield.

SPECIFICATIONS

Model	BKB-5601
Max. Speed	4000rpm
Max. RCF	3150xg
Max. Capacity	4x100ml
Speed Accuracy	±10r/min
Time Setting Range	1min to 99min
Temperature Setting Range	Room temperature+10°C to 70°C
Noise	≤65dB(A)
Power Supply	AC220V±22V 50Hz
Power Consumption	1500W
Dimension (WxDxH)	520x550x400mm
Packing Dimension (WxDxH)	772x772x582mm
Weight	85kg
Alt Name	Oil Testing Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Swing Rotor Yoke	4000	3150	36x10ml
2	Swing Buckets	4000	3150	36x10ml
3	Swing Rotor Yoke	1700	691	4x100ml
4	Swing Buckets	1700	691	4x100ml

FEATURES

LCD of centrifuge machine for crude oil testing to show rotor No., speed, time, temperature, RCF, Acc/Dec rate, and trouble information/error code.

A circular heater for fast heating and temperature uniformity realizes the functions of temperature controlling and temperature holding, which meet the temperature requirements for the oil industry, making it a perfect centrifuge machine for crude oil testing. The frequency-conversion motor of centrifuge machine for crude oil testing with high torque is characterized by stable and fast acceleration/deceleration as well as long service life, maintenance-free, and low noise for stable operation.

The centrifuge for crude oil testing is equipped with a specially designed tube, making it resistant to damage during operation. The oil centrifuge is a widely utilized tool in Daqing Oilfield.

The crude oil centrifuge for sale at Cence features a self-diagnosis system that automatically detects faults such as over-speed and over-temperature, displaying error codes on the screen to enhance safety performance.

The centrifuge machine for crude oil testing offers simple operation with microprocessor control and a user-friendly touch panel.

Parameters like speed, RCF, temperature, time, and Acc/Dec rates are fully adjustable. The parameters of the centrifuge machine for crude oil testing can be revised during a run for convenient use.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com