



OIL TEST CENTRIFUGE BKB-5701

OIL TEST CENTRIFUGE BKB-5701

OIL TESTING CENTRIFUGE



LCD to show rotor No., speed, time, temperature, RCF, ACC/DEC rate, and error information/error code.

With the function of temperature-control and constant, the oil centrifuge adopts a band heater, fast heating speed, and temperature uniform.

Stainless steel chamber, safe and reliable.

The glass oil tube is manufactured by a special process, which won't be broken so easily during operation, and it's widely used in Daqing Oilfield.

SPECIFICATIONS

Model	BKB-5701
Max. Speed	4000rpm
Max. RCF	3150xg
Max. Capacity	4x100ml
Speed Accuracy	±10r/min
Time Setting Range	1min to 99min
Temperature Setting Range	Room temperature+10°C to 70°C
Noise	≤65dB(A)
Power Supply	AC220V±22V 50Hz
Power Consumption	1500W
Dimension (WxDxH)	560x660x850mm
Packing Dimension (WxDxH)	756x856x1320mm
Weight	108kg
Alt Name	Oil Testing Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Swing Rotor Yoke	4000	3150	36x10ml
2	Swing Buckets	4000	3150	36x10ml
3	Swing Rotor Yoke	1700	691	4x100ml
4	Swing Buckets	1700	691	4x100ml



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com