



BLOOD BANK CENTRIFUGE BKB-5901

BLOOD BANK REFRIGERATED CENTRIFUGE BKB-5901

LOW SPEED BLOOD BANK CENTRIFUGE



Adopting microcomputer control, LCD display, brushless DC motor drive, and 10 user programs for storage.

Short acceleration/deceleration time, with a minimum acceleration time of $\leq 15S$. Equipped with a fault self-diagnosis system, this blood centrifuge machine can automatically detect various faults such as overspeed, imbalance, and door cover self-locking, and display fault information in the display window to improve the safety performance of the instrument.

Can store up to 50 user-defined programs.

Equipped with inching function, it can quickly complete separation.

Electronic door lock, stainless steel inner chamber, humanized design, with minimal force when closing the door cover. Gently press the door cover to automatically lock.

SPECIFICATIONS

Model	BKB-5901
Max. Speed	5000rpm
Max. RCF	3555xg
Max. Capacity	12x15ml
Speed Accuracy	$\pm 30r/min$
Time Setting Range	1min to 99min
Noise	AC220V $\pm 22V$ 50/60Hz
Power Supply	$\leq 65dB(A)$
Power Consumption	250W
Dimension (WxDxH)	380x500x300mm
Packing Dimension (WxDxH)	465x660x350mm
Weight	19kg
Alt Name	Low Speed Blood Bank Centrifuge

ITEMS INCLUDED

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	5000	3555	12*15 ml



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com