



CELL WASHING CENTRIFUGE BKB-6001

CELL WASHING CENTRIFUGE BKB-6001

CELL WASHING CENTRIFUGE



Special damping structure to reduce vibration.
Integrated molding of housing.
Microprocessor control, DC brush-less motor.
Touch panel and LCD display.
Countdown time displayed by seconds if the centrifugation time is less than 1min.
Electronic door lock, centrifuge cannot work if the door is open and the door cannot be opened when it is working.

SPECIFICATIONS

Model	BKB-6001
Max. Speed	4500rpm
Max. RCF	2000xg
Max. Capacity	(0.25~1ml)x12
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min
Noise	≤60dB(A)
Power Supply	AC220V ±22V 50/60Hz
Power Consumption	100W
Dimension (WxDxH)	456x330x266mm
Packing Dimension (WxDxH)	568x423x300mm
Weight	12kg
Alt Name	Cell Washing Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	SERO Rotor	3000	1000	12x7ml
2	HLA Rotor	4500	2000	0.25~1ml glass tube 12pcs



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com