



## HEMATOCRIT CENTRIFUGE BKB-6101

# HEMATOCRIT CENTRIFUGE BKB-6101

## HEMATOCRIT CENTRIFUGE



Parameters can be set freely during a run.

Designed for clinical efficiency, our hematocrit centrifuge machine incorporates microprocessor control and a DC brushless motor for stable and low noise operation in busy medical environments.

G-force can be calculated automatically based on the set hematocrit centrifuge rpm, simplifying the operation for users.

With a self-diagnosis system, this blood centrifuge machine can automatically identify faults, including overspeed, no speed signal, and door lid interlock issues, significantly improving safety.

The countdown time is displayed in seconds if the centrifugation time is less than 1 minute.

5-inch HD touch screen, simple interface, setting parameters, and running parameters at a glance.

## SPECIFICATIONS

|                           |                       |
|---------------------------|-----------------------|
| Model                     | BKB-6101              |
| Max. Speed                | 12000rpm              |
| Max. RCF                  | 14800xg               |
| Max. Capacity             | 36 capillaries        |
| Speed Accuracy            | ±30/min               |
| Time Setting Range        | 1min to 99min         |
| Noise                     | ≤ 65dB(A)             |
| Power Supply              | AC220V±22V 50Hz       |
| Power Consumption         | 300W                  |
| Dimension (WxDxH)         | 320x370x235mm         |
| Packing Dimension (WxDxH) | 530x410x290mm         |
| Weight                    | 11kg                  |
| Alt Name                  | Hematocrit Centrifuge |

## ACCESSORIES FOR PURCHASE

| No | Name             | Max. Speed | Max. RCF | Max. Capacity  |
|----|------------------|------------|----------|----------------|
| 1  | Hematocrit rotor | 12000      | 14800    | 12 capillaries |
| 2  | Hematocrit rotor | 12000      | 14800    | 24 capillaries |
| 3  | Hematocrit rotor | 12000      | 14800    | 36 capillaries |

## FEATURES

Parameters can be set freely during a run.

Designed for clinical efficiency, our hematocrit centrifuge machine incorporates microprocessor control and a DC brushless motor for stable and low noise operation in busy medical environments.

G-force can be calculated automatically based on the set hematocrit centrifuge rpm, simplifying the operation for users.

With a self-diagnosis system, this blood centrifuge machine can automatically identify faults, including overspeed, no speed signal, and door lid interlock issues, significantly improving safety.

The countdown time is displayed in seconds if the centrifugation time is less than 1 minute.

5-inch HD touch screen, simple interface, setting parameters, and running parameters at a glance.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)